

UpDate!



1896 SCHWERT
SWORD SYNONYM FÜR QUALITÄT & VERTRAUEN
 SYNONYM FOR QUALITY & TRUST



Weingart Extended Ligaturenzange Weingart Extended Ligature Pliers

Modern Therapy · Modern Instruments

SCHWERT
 SWORD · ESPADA
 EPEE · SPADA



5308-L



SCHWERT
SWORD · ESPADA
EPEE · SPADA

UpDate!

Weingart Extended Ligaturenzange · Weingart Extending Ligature Pliers

Modern Therapy · Modern Instruments

Mit der neuen SCHWERT Extended Weingart Ligaturenzange, mit einer Länge von 155 mm, ist der Forderung nach einem besonders flexibel einsetzbaren Instrument für die KFO Rechnung getragen worden. Die längeren Griffe erlauben einen optimalen Zugang zu schwer erreichbaren Regionen, wie es beispielsweise in der Lingual-technik öfter der Fall ist. Sie gewährleisten außerdem ein taktiles und ausgewogenes Arbeiten. Die grazile Ausarbeitung und präziser Grip des Arbeitsendes in Kombination mit einem hochwertigen Oberflächenfinish sorgen für maximalen Komfort für Anwender und Patienten.

Die SCHWERT Zange ist mit einem gerieften Hartstahleinsatz (HCS) versehen, welcher eine lange Lebensdauer des Instruments und die bestmögliche Adaption in den unterschiedlichen Anwendungsbereichen sicherstellt.

Durch die bessere Erreichbarkeit von schwer zugänglichen Bereichen, können Behandlungszeiten verkürzt und somit höhere Wirtschaftlichkeit in der Praxis erreicht werden.



5308-L Weingart Extended

5309 Weingart Standard

With the new SCHWERT Extended Weingart ligature pliers, with an extended length of 155 mm, an instrument very flexible in its use for orthodontics was created. The longer handles allow an optimal approach especially in hardly reachable areas, e.g. in the lingual-technique. Furthermore they allow a very tactile and balanced work. Due to the fine working tip and its perfect grip, in combination with a high-quality finish, the pliers provide maximum comfort for both, patient and operator.

The pliers tips are made of high carbon steel (HCS) inserts, which extends the durability considerably and ensures the best adaption possible in every area.

Due to the easier reachability of hard approachable areas, the instrument helps to shorten treatment times and supports a better profitability in the practice.

HCS

Hartstahl Einsätze zeichnen sich durch ihre Härte und Langlebigkeit aus.

High Carbon Steel Inserts are known for their hardness and durability.



5308-L

15,5 cm, • max 0,7 mm / .028", HCS-Einsatz
15.5 cm, • max 0.7 mm / .028", HCS-insert